

EV Analysis Platform

Nanoparticle Tracking Analysis Report

Video Analysis | Laser Scattering Microscopy

Operator: System | Report Generated: 2026-01-20 12:38

Sample Parameters


Sample Name:	PC3_100kDa_F9T15	Result (sizes in nm)		
Sample Info:	PC3 100kDa Fraction - 20251217_0019	Number	Volume	
Electrolyte:	WATER	Median (X50)	162.5	162.5
Temperature:	25.0 °C	StdDev	83.1	83.1
Concentration:	7.60E+07 Particles/mL			
Dilution Factor:	500			

Instrument Parameters

Laser Wavelength:	488 nm	Frame Rate:	30 fps
Filter:	Scatter	Video Resolution:	medium
Sensitivity:	80	Traced Particles:	152
Shutter:	100	Positions:	11

Size Distribution

Size Distribution - PC3_100kDa_F9T15



X Values (all sizes given in nm)

	Number	Concentration	Volume
X10	107.5	107.5	107.5
X50	162.5	162.5	162.5
X90	297.5	297.5	297.5

Span	1.2	1.2	1.2
Mean	182.7	182.7	182.7
StdDev	83.1	83.1	83.1

Peak Analysis

Diameter (nm)	Particles/mL	FWHM (nm)	Percentage
122.5	7.60E+07	83.1	100.0

Analysis Software: EV Analysis Platform v1.0 | Report Generated: 2026-01-20 12:38:30

Comments:

_____ (Signature)